Package insert Testsealabs Covid-19 Antigen (SARS-CoV-2) Test Cassette (Swab)

Rapid test for qualitative! Detection of SARS-CoV-2 nucleocapsid antigen in nasopharyngeal, oropharyngeal and nasal swabs. For professional in vitro diagnostic use only

INTENDED USE

The COVID-19 Antigen Test Cassette is a chromatographic immunoassay for the qualitative detection of SARS-CoV-2 nucleocapsid antigen in nasopharyngeal, oropharyngeal and nasal swabs to aid in the diagnosis of SARS-CoV-2 virus infection.

PRINCIPLE

The COVID-19 Antigen Test Cassette is a qualitative membrane strip based immunoassay for the detection of COVID-19 antigen in nasopharyngeal, oropharyngeal and nasal swab specimen. In this test procedure, anti-SARS-CoV-2-N antibody is immobilized in the test line region of the device. After a nasopharyngeal swab specimen is placed in the specimen well, it reacts with anti-SARS-CoV-2-N antibody coated particles that have been applied to the specimen pad. This mixture migrates chromatographically along the length of the test strip and interacts with the immobilized anti-SARS-CoV-2-N antibody. If the specimen contains SARS-CoV-2 antigen, a colored line will appear in the test line region indicating a positive result. If the specimen does not contain SARS-CoV-2 antigen, a colored line will not appear in this region indicating a negative result. To serve as a procedural control, a colored line will always appear at the control line region indicating that proper volume of specimen has been added and membrane wicking has occurred.

REAGENTS

The test contains anti-SARS-CoV-2-N antibody as the capture reagent, another anti-SARS-CoV-2-N antibody as the detection reagent. A Goat anti-Mouse IgG is employed in the control line system.

PRECAUTIONS

- For professional in vitro diagnostic use only. Do not use after the expiry date.
- Make sure that the foil pouch containing the test cassette is not damaged before opening it for use.
- Wear gloves and personal protective equipment when taking and applying samples. Do not touch the reagent membrane and the sample window.
- Do not eat, drink or smoke in the area where specimens and kits are handled.
- Handle all samples as if they contained infectious agents.
- Observe established precautions against microbiological hazards throughout all procedures and follow the standard procedures for proper disposal of specimens.
- Conduct the test at a room temperature of 15 - 30 °C. Humidity and temperature can adversely affect the results.

STORAGE AND STABILITY

Store as packaged in the sealed pouch at room temperature or refrigerated (4-30°C).

The test is stable to the expiration date printed on the sealed pouch. The test must remain in the sealed pouch until use. DO NOT FREEZE. Do not use beyond the expiration date.

SPECIMEN COLLECTION AND PREPARATION

The COVID-19 Antigen Test Cassette is designed for use with nasopharyngeal, oropharyngeal and nasal swabs. Have the swab performed by a medically trained person. For best results, a nasopharyngeal swab is recommended.

Instructions for the nasopharynx swab procedure

Insert the swab liberally through the nose into the nasopharynx and swab the nasopharynx in 2-3 circular motions.

Instructions for the nasal swab procedure

Insert the entire tip of the swab two to three centimetres into the left nostril. Rub the inside of the nostril in a circular motion for at least 15 seconds. Remove the swab and insert it into the right nostril. Swab the inside of the nostril in a circular motion for at least 15 seconds.

General information

Do not return the swab to its original paper wrapper. For best results, swabs should be tested immediately after collection. If it is not possible to test immediately, it is strongly recommended that the swab is placed in a clean, unused plastic tube labelled with patient information to maintain best performance and avoid possible contamination. The sample can be kept tightly sealed in this tube at room temperature (15-30°C) for a maximum of one hour. Make sure that the swab is firmly seated in the tube and that the cap is tightly closed. If a delay of more than one hour occurs, discard the sample. A new sample must be taken for the test.

If specimens are to be transported, they should be packaged according to local regulations for the transport of atiological agents.

MATERIAL

Materials provided

Test device Extraction tubes with sample buffer
Package insert Sterile swab Workstation

Materials required but not provided: Timer

DIRECTIONS FOR USE

Allow the test, sample and buffer to reach room temperature 15-30°C (59-86°F) before running.

- Place the extraction tube in the workstation.
- Open the extraction tube filled with sample buffer and remove the cap
- Have the nasopharyngeal, oropharyngeal or nasal swab carried out by a medically trained person as described.
- Place the swab in the extraction tube. Rotate the swab for about 10 seconds while dripping the head against the inside of the tube to release the antigen in the swab.
- Remove the swab while pressing the head of the swab against the inside of the extraction tube to expel as much liquid as possible from the swab. Dispose of the swab according to the disposal regulations for biological waste.
- Twist the cap on the extraction tube and make sure it is tight. Now carefully twist off the cap over the filter valve.
- Place 3 drops of the sample vertically into the sample window of the test cassette. Read the result after 10-15 minutes. Read the result within 20 minutes. Otherwise, a repetition of the test is recommended.

EVALUATION OF RESULTS


Positive: Two red lines appear. One red line appears in the control zone (C) and one red line in the test zone (T). The test is considered positive if even a faint line appears. The intensity of the test line can vary depending on the concentration of the substances present in the sample.

Negative: Only in the control zone (C) a red line appears, in the test zone (T) no line appears. The negative result indicates that there are no SARS-CoV-2 antigens in the sample or the concentration of the antigens is below the detection limit.

Invalid: No red line appears in the control zone (C). The test is invalid even if there is a line in the test zone (T). Insufficient sample volume or incorrect handling are the most likely reasons for failure. Review the test procedure and repeat the test with a new test cassette.

QUALITY CONTROL

The test contains a coloured line that appears in the control zone (C) as an internal

procedural control. It confirms sufficient sample volume and correct handling. Control standards are not supplied with this kit. However, it is recommended that positive and negative controls be tested as good laboratory practice to confirm the test procedure and verify proper test performance.

LIMITATION

- This test detects both viable (live) and non-viable, SARS-CoV and COVID-19. Test performance depends on the amount of virus (antigen) in the sample and may or may not correlated with viral culture results performed on the same sample.
- A negative test result may occur if the antigen concentration in a sample is below the detection limit of the test. The detection limit of the test was determined with recombinant SARS-CoV-2 nucleocapsid protein and is 100 pg/ml.
- The performance of the SARS-CoV-2 Antigen Test Cassette has been evaluated using only the procedures described in this package insert. Modifications to these procedures may alter the performance of the test.
- False negative results may occur if a sample is collected, transported or handled improperly.
- False results may occur if samples are tested later than one hour after collection. Samples should be tested as soon as possible after collection.
- Positive test results do not exclude co-infection with other pathogens.
- Negative test results are not intended to be predictive of viral or bacterial infections other than SARS-CoV-2.
- Negative results from patients with symptom onset after more than seven days should be treated as presumptive and confirmation with another molecular assay should be performed.
- When differentiation of specific SARS-CoV-2 strains is required, additional testing is required in consultation with state or local health authorities.
- Children may tend to shed viruses longer than adults, which can lead to differences in susceptibility between adults and children and more difficult comparisons.

PERFORMANCE CHARACTERISTICS

Detection limit: The detection limit of the test was determined with infectious SARS-CoV-2 virus and is 50 x TCID₅₀.

Clinical Performance Characteristics of Nasopharyngeal Swab

To determine the sensitivity and specificity, the COVID-19 Antigen Test Cassette was compared with nasopharyngeal swabs using a commercial PCR test.

Sensitivity	97.6% 95% CI: (94.9%-100%)
Specificity	98.4% 95% CI: (96.9%-99.9%)

Determination of the specificity

Number of samples	PCR negative	COVID-19 Antigen Test Cassette
250	250	246/250 = 98.4%
Total	250	246/250 = 98.4% 95% CI: (96.9%-99.9%)

Determination of the sensitivity

Days post Symptom onset	Number of samples	PCR positive	COVID-19 Antigen Test Cassette
1	3	3	3/3=100%
2	9	9	9/9=100%
3	14	14	14/14=100%
4	17	17	17/17=100%
5	22	22	21/22=95.4%
6	26	26	25/26=96.1%
7	34	34	33/34= 97.0%
Total	125	125	122/125=97.6% 95% CI: (94.9%-100%)

Clinical Performance Characteristics of Nasal swab

To determine the sensitivity and specificity, the COVID-19 antigen rapid test kit with nasal swabs was compared with a commercial PCR test. All swabs were taken by the patients themselves and the results were evaluated by the patients together with the nursing staff.


Sensitivity	93.6% 95% CI: (92.5%-94.7%)
Specificity	98.8% 95% CI: (98.5%-99.1%)

Determination of the specificity

Number of samples	PCR negative	COVID-19 Antigen Test Cassette
250	250	247/250 = 98.8%
Total	250	247/250 = 98.8% 95% CI: (98.5%-99.1%)

Determination of the sensitivity

Days post Symptom onset	Number of samples	PCR positive	COVID-19 Antigen Test Cassette
1	3	3	3/3=100%
2	9	9	9/9=100%
3	14	14	14/14=100%
4	17	17	17/17=100%
5	22	21	21/22=95.4%
6	26	26	24/26=92.3%
7	34	34	29 /34= 87.8%
Total	125	125	117/125=93.6% 95% CI: (92.5%-94.7%)

Cross-reactivity: The COVID-19 antigen rapid test cassette has been tested for specificity and cross-reactivity with other pathogens that may cause similar symptoms. The results showed no cross-reactivity.

Pathogen	Concentration
Pseudomonas aeruginosa	1x10 ⁸ org/mL
Streptococcus sp group F	1x10 ⁸ org/mL
Streptococcus salivarius	1x10 ⁸ org/mL
Streptococcus pyogenes	1x10 ⁸ org/mL
Streptococcus pneumoniae	1x10 ⁸ org/mL
Staphylococcus epidermidis	1x10 ⁸ org/mL
Staphylococcus aur. subsp. aureus	1x10 ⁸ org/mL
Nisseria subflava	1x10 ⁸ org/mL
Neisseria lactamica	1x10 ⁸ org/mL
Moraxella catarrhalis	1x10 ⁸ org/mL
Escherichia coli	1x10 ⁸ org/mL
Corynebacterium	1x10 ⁸ org/mL
Candida albicans	1x10 ⁸ org/mL
Arcanobacterium	1x10 ⁸ org/mL
Human Coronavirus OC43	2.45x 10 ⁶ LD ₅₀ /ml
Human Coronavirus NL63	1.17x10 ⁶ U/ml
Influenza A H1N1	3.16 x 10 ⁶ TCID ₅₀ /ml
Influenza A H3N2	1 x 10 ⁶ TCID ₅₀ /ml
Influenza B	3.16x10 ⁶ TCID ₅₀ /ml
Human Rhinovirus 2	2.81 x10 ⁶ TCID ₅₀ /ml
Human Rhinovirus 14	1.56x 10 ⁶ TCID ₅₀ /ml
Human Rhinovirus 16	8.89x10 ⁶ TCID ₅₀ /ml
Masem	1.56x10 ⁶ TCID ₅₀ /ml
Mumps	1.56x 10 ⁶ TCID ₅₀ /ml
Parainfluenza Virus 2	1.56x 10 ⁷ TCID ₅₀ /ml
Parainfluenza Virus 3	1.58 x 10 ⁸ TCID ₅₀ /ml
Respiratorisches Syncytial-Virus	8.89x10 ⁴ TCID ₅₀ /ml

Interfering substances: The following compounds were tested with the COVID-19 rapid antigen test and no interferences was observed.

Substance	Concentration	Substance	Concentration
whole blood	20µl/ml	Mupirocin	12 mg/ml
Muzin	50µg/ml	Oxymetazoline	0.6 mg/ml
Budesonid Nasenspray	200µl/ml	Phenylephrine	12 mg/ml
Dexamethasone	0.8 mg/ml	Rebetol	4.5µg/ml
Flunisolide	6.8 ng/ml	Relenza	282 ng/ml

Symbol	Meaning	Symbol	Meaning
	Medical in vitro diagnosis		Storage temperature limits (4-30°C)
	Manufacturer		Tests per set
	Batch code		Do not reuse
	Follow the instructions for use		Authorised Representative in the European Community
	Expiry date		Catalogue number

HANGZHOU TESTSEA BIOTECHNOLOGY CO.,LTD.
3rd Floor, Building 6,
No.8-2 Keji Road, Yuhang District,
Hangzhou, China

EC REP Lotus NL B.V.
Koningin Julianaplein 10,
1e Verd, 2595AA, The Hague,
Netherlands.

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